


Remove fixation. Reveal the nystagmus.

+20 D binocular, rechargeable, bedside-ready Frenzel goggles.



WHY EQUIGLASSES



+20 D binocular

High-diopter lenses magnify the eyes while blurring patient vision.



Built-in illumination

Eyes stay brightly visible even in a darkened room.



Rechargeable

Cordless — no disposable cells, no trailing wires.



Bedside-ready

Lightweight, point-of-care form for clinic and ER.

FIVE BEDSIDE EXAMS. ONE DEVICE.

1

Spontaneous Nystagmus

Resting nystagmus with fixation abolished.

2

Gaze-Evoked Nystagmus

Eccentric gaze under fixation removal.

3

Positional Testing

Dix-Hallpike, Roll Test and bedside manoeuvres.

4

Head-Shake Test

Post-headshake nystagmus with fixation suppressed.

5

Provocative Manoeuvres

Valsalva, hyperventilation, vibration and Tullio.